TECHNICAL CHARACTERISTICS

- Response time: : >1 to <5 min
- Working Temperature: 5-50°C
- Minimum required volume : 0.5 mL

TISAB* (10%): For samples with high ionic strength (>0.01M), pH out of optimum range or with presence of complexing ions (Al³⁺, Cu²⁺, Fe³⁺), use the special buffer: 17.2g Sodium Acetate (anhydrous), 8g Sodium Chloride, 1g CDTA and 2.7ml Acetic Acid (98%) in 250ml of DI water.

**10% of TISAB is usually enough for buffering samples*

*Supporting information about Fluoride measurement is online available:
http://www.ntsensors.com/combis/C19_support.pdf*

ADDITIONAL EQUIPMENT

- NT ION METER or an equivalent meter: pH/mV-meter with resolution of 0.1mV.
- Flasks and pipettes.

REAGENTS

- Deionized water.
- Standard and conditioning solutions of the primary ion to be determined.

PREPARATION AND USE

Before using the CNT ISE COMBI, it is recommended to read the instructions of your meter.

Condition the CNT ISE in a solution of the target ion of 1000ppm at least for 10 minutes before daily use.

If the sensor is new, has been prolonged time without use, or has been in contact with interference containing samples, condition time is

recommended to be minimum 8 hours or until stable potential reading.

If target concentration is lower than 100ppm, is recommended a second conditioning process in 100ppm (even less) solution at least for 10 minutes.

1. Plug the BNC terminal of the CNT ISE to the meter.
2. Remove the protection cap from the sensing area.
3. Calibrate the sensor ⁽¹⁾.

(1)Regarding the complexity of the sample matrix and some different factors, the analytical procedure could be direct calibration or different analytical techniques, such as the standard addition, etc.

4. Rinse with DI water and dry the outer body with a clean tissue.
5. Measure the sample.
6. Rinse with DI water and dry the outer body of the probe between each sample measure.
7. Keep dry and clean with the protective cap.

Presence of solid particles in suspension and turbid solutions do not affect to the overall performance of the sensor.

RECOMENDATIONS

- ✓ Keep constant the same conditions of temperature, stirring, both in samples and standards.
- ✓ Follow the instructions for better conservation of the sensor.
- ✓ Great care has to be taken to do not damage the tip. The electrode can be irreversibly damaged if this part is hit or grated.

GUARANTY

Sensors are guaranteed of any manufacturing defect.

NT Sensors will replace without additional cost the Sensors which, after being revised by its technical post-service have been considered as "defect from origin".


The Guaranty of the sensors does not cover the defects caused by:

- inadequate use,
 - the usual aging of the sensor,
 - the logic premature aging caused by certain samples,
 - the damaged caused by accident.
- The guaranty is valid through a period of 6 months.

MAINTENANCE AND STORAGE

- ✓ The CNT ISE COMBI does not require maintenance due to not contain internal liquid solutions.
- ✓ Place the protective cap when not use the sensor. Do not leave the sensing area in contact with air/atmosphere for longer time than necessary.
- ✓ Storage at temperatures below 50°C.
- ✓ Storage in a dry, cool place avoiding the direct contact with the sunlight.

Specifications Especificaciones	F ⁻
Slope (mV/dec) at 25°C <i>Pendiente a 25°C</i>	-54 ± 5
Range (mg/L) <i>Rango</i>	0.1 to 1900
Range (mol/L) <i>Rango</i>	5x10 ⁻⁶ to 0.1
pH Range <i>Rango pH óptimo</i>	4 to 8
Main interferences (log K _{ij}) <i>Interferentes principales</i>	OH ⁻ (-1)* *Kept pH below 8 *Mantener pH por debajo de 8

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http://www.ntsensors.com/combi/C19_support.pdf

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